



**CAPT\_QUA\_0700 Glucose Meter Proficiency Testing**

|             |               |               |                 |
|-------------|---------------|---------------|-----------------|
| Identifier: | CAPT_QUA_0700 | Version #:    | 1.6             |
| Folder:     | POCT W.BLOOD  | Type:         | Procedure (2yr) |
| Subfolder:  | QC & PT       | Effective on: | 2025-02-24      |

# Glucose Meter Proficiency Testing

# CAPT\_QUA\_0700

**Purpose:**

- The purpose of this document is how to perform proficiency testing on the Nova glucose meters.
- All POCT is evaluated using a formal proficiency testing program or by an alternate means of assessment.
- POCT PT assessment is examined by personnel who routinely examine patient samples. Proficiency testing or alternate assessment is performed by POCT operators.
- Proficiency testing or alternate assessment occurs at least twice per year.

*PT Frequency*

*DAP Provisionally Accredited Facility*

|                           |              |   |
|---------------------------|--------------|---|
| DAP reportable measurands | All services | Minimum two samples and one test event prior to full DAP assessment |
| Non-reportable measurands |              |   |

*DAP Accredited Facility*

|                           |              |  |
|---------------------------|--------------|--|
| DAP reportable measurands | All services | Minimum four samples per year<br>Minimum two testing events per year |
| Non-reportable measurands |              |  |

- Proficiency control material for quality assurance purposes.
- The alternative assessment is established by the laboratory medical director or designate.
- Proficiency testing or alternate assessment results are monitored by the laboratory medical director or designate at a defined interval and discussed with relevant personnel.
  - Designate reviews all PT results.
  - If non-conforming results that require investigation, will be brought to attention of medical director.
- The laboratory does not communicate with other participants in the inter-laboratory comparison (PT) program about sample data until after the date for submission of the data. The laboratory does not refer inter-laboratory (PT) comparison samples for confirmatory examinations before submission of the data, although this would routinely be done with patient samples.

|              |               |                      |               |
|--------------|---------------|----------------------|---------------|
| Written by:  | Ron Garbuio   | Approved by (sign.): |               |
| Reviewed by: | Ron Garbuio   | Approved by (name):  | Cheng-Han Lee |
| Reviewed on: | 2019-08-27    | Approved on:         | 2019-08-27    |
| Renewed by:  | Ronny Garbuio | Revision Date:       | 2027-02-24    |
| Renewed on:  | 2025-02-24    |                      |               |

**Procedure:****1. Ordering Proficiency Material:**

Obtain proficiency material from reputable source. Obtained from College of American Pathologists.

**QUALITY CROSS CHECK-WHOLE BLOOD GLUCOSE - WBGQ**

Analytes/procedures in **bold** type are regulated for proficiency testing by the Centers for Medicare & Medicaid Services (CMS).

| Analyte | Challenges/Shipment | Number of Shipments    |
|---------|---------------------|------------------------|
| Glucose | 3                   | Two shipments per year |

The CAP Accreditation Program requires all accredited laboratories performing waived whole blood glucose testing using glucose meters to perform alternative performance assessment. This program can be used to meet alternative performance assessment requirements.

**Additional Information**

**Quality Cross Check** is a convenient solution to monitor instrument performance and assess comparability across multiple instruments in your laboratory and to identify potential issues before they affect patient results.

**How It Works:**

Receive three challenges in each of two mailings a year.  
Report up to 30 instruments for each challenge.  
Receive a custom report package that includes peer group comparison and instrument comparability statistics for each reported analyte.

**Program Information**

Three 2.0-mL whole blood specimens  
Report up to 30 instruments

**Shipping Schedule**

Shipment A: April  
Shipment B: September

## 2. Prepare Proficiency Material:

### a. Supplies:

#### i. Proficiency Material:

1. 3 vials
2. Store at room temperature; do not freeze
3. Viable one month after opening

#### ii. Micro-Cuvettes: 3 per glucometer (45 total)

#### iii. Micro-Cuvette holder: borrow from cancer genetics

#### iv. Micro-Pipettes: 1 per cuvette (45 total)

#### v. Biohazard bag: 1 per glucose meter (15 total)

#### vi. Labels: 1 per cuvette (see part b.)

#### vii. Bi-Yearly Proficiency Report Template: 1 per meter

Refer to: [CAPT\\_QUA\\_0700c11 Bi-Yearly Proficiency Report Template](#)

### b. Label each cuvettes with proficiency test codes ex. 01, 02, 03 etc

### c. In biological safety cabinet:

#### i. Wear personal protective equipment

1. Gown
2. Gloves

#### ii. Perform transfer one PT sample at a time

#### iii. Place (15) cuvettes in tube holder

#### iv. Add minimum 2 - 3 drops of PT per cuvette

#### v. Cap.

#### vi. Repeat for PT 2 and PT 3

#### vii. Place one of each PT sample cuvette into a biohazard bag with 3 micropipettes. Place template proficiency report in sleeve of biohazard bag.

#### viii. Remove personal protective equipment

### d. Delivery proficiency testing to wards/locations.

**3. Analyze this specimen as you would a patient's blood sample. Instructions:**

- Wear gloves
- Check Meter Name/Serial Number
- Prepare Nova StatStrip Glucose Meter for Patient Testing. \*Perform QC if not already performed within 24hrs.
- **Gently mix the sample.** Tap sample to the bottom of the test tube.
- Apply sample to the end of the test strip using plastic dropper provided.
- Record the Glucose Meter results in the area provided on the template.
- Record test strip code and date & time of testing. Sign initials & operator ID.
- Dispose of samples in garbage or biohazard container.
- Send results to:
  1. Directly Labs Site Supervisor (room 3208A)
  2. Scan and email ([rgarbuio@bccancer.bc.ca](mailto:rgarbuio@bccancer.bc.ca)) with attention to Ron Garbuio, POCT PT Results.
  3. Fax. Attn: Ron Garbuio 604 877-6178

**4. Send the results to CAP online. Do not fax or email. – done by POCT coordinator**

- a. Log into the CAP website
- b. Select View, Enter or submit PT results

**Welcome to the CAP**

View, enter, or submit PT Results →

---

Claim CME/CE credit for faxed AP results →

---

Access the Cancer Protocols →

- c. Find the proficiency testing in the list

| Due Date   | Ship Date  | Mailing     | Kit #    | Seq | Status   | CAP #      | Owned by |
|------------|------------|-------------|----------|-----|----------|------------|----------|
| 04/16/2019 | 03/25/2019 | CYP-A 2019  | 31966702 | 1   | Received | 7190232-01 |          |
| 04/23/2019 | 04/01/2019 | WBGQ-A 2019 | 32264243 | 1   | Received | 7190232-01 |          |

**Step 1: Enter/Edit/View results**

Total of 6 page(s)

| Page | Status       | Date Received    | Via    | Data                       |
|------|--------------|------------------|--------|----------------------------|
| 1    | Received     | 4/11/19 12:20 PM | Online | <a href="#">View/Edit</a>  |
| 2    | Received     | 4/11/19 12:20 PM | Online | <a href="#">View/Edit</a>  |
| 3    | Received     | 4/11/19 12:20 PM | Online | <a href="#">View/Edit</a>  |
| 4    | Not Received |                  |        | <a href="#">Enter Data</a> |
| 5    | Not Received |                  |        | <a href="#">Enter Data</a> |

British Columbia Cancer Agency  
Vancouver, BC V5Z 4E6

**OTHER ACTIONS**

- [Print a blank result form](#)
- [View/print your saved data](#)
- [View kit transaction history](#)
- [View kit instructions](#)
- [Contact Us](#)

**Step 2: Review and submit results**

Results have been submitted and received successfully.

APPROVE AND SUBMIT TO CAP

- d. Enter all the results:
  - i. Unit of measurement: mmol/L
  - ii. Enter the site

|   |                        |
|---|------------------------|
| A | 6E CC                  |
| B | 6E PT                  |
| C | 5E CC                  |
| D | 5E PT                  |
| E | 5W PT                  |
| F | 3 LABS-DI AMB          |
| G | 2 AMB                  |
| H | FAIRMONT               |
| I | 2 RADIATION THERAPY    |
| J | 2 RADIATION THERAPY CC |
| K | SURGICAL DAYCARE       |
| L | PET1                   |
| M | PET2                   |
| N | PET3                   |

- iii. Exceptions – with problem PT
  - iv. Results for each PT
  - v. Save and submit
- e. Save copy in: [H:\Lab\\_Medicine\Pathology Office\Site Supervisor\POCT\Proficiency Testing - CAP\Results](H:\Lab_Medicine\Pathology Office\Site Supervisor\POCT\Proficiency Testing - CAP\Results)
- f. Enter results in spreadsheet: [H:\Lab\\_Medicine\Pathology Office\Site Supervisor\POCT\Proficiency Testing – CAP\Results\Proficiency Testing Results.xlsx](H:\Lab_Medicine\Pathology Office\Site Supervisor\POCT\Proficiency Testing – CAP\Results\Proficiency Testing Results.xlsx)


## 5. POCT PT Assessments

- a. Results are monitored by the laboratory medical director or designate within 4 weeks of receiving and discussed with relevant personnel; ex. Point of Care Committee Meetings and Quality Improvement Monthly Meetings.

Proficiency testing or alternate assessment results are monitored by the laboratory medical director or designate at a defined interval and discussed with relevant personnel.


- i. If all acceptable with no impact to patient care, Site Supervisor/POCT Quality Coordinator can sign the final report

- ii. If non-conforming results, investigated and have signed by Medical Director
- b. Obtain report
- i. Go to the CAP website: <https://www.cap.org>
  - ii. Login
  - iii. Got to Access e-LAB Solutions Suite
  - iv. Select Evaluation Reports



## Welcome to the CAP

- v. Select the proficiency test, select the report, and view report



- vi. Save a copy of the report in: [H:\Lab\\_Medicine\POCT - Point of Care\Proficiency Testing - CAP\Results](H:\Lab_Medicine\POCT - Point of Care\Proficiency Testing - CAP\Results)
- c. Unacceptable POCT PT assessment results are investigated and corrective action is implemented where indicated. This investigation and any corrective action is documented and retained.
- i. Laboratory will try to repeat test with same sample on same meter to determine if the error was operator or meter related.
  - ii. Implement corrective actions

- iii. Check for trends that indicate potential nonconformities
  - iv. Analyse the root cause – check the PT assessment review/recommendation and consult pathologists
  - v. Take further action as required to prevent occurrence
  - vi. File the record of all corrective action
- d. The authority to withdraw equipment or discontinue a POCT examination in the event of serious POCT PT or alternate assessment problems is defined.
- i. a clinically significant impact to patient results has been confirmed
  - ii. the accuracy and reliability of test results cannot be verified, or
  - iii. the cause of significant or ongoing PT exceptions cannot be determined
- e. When predetermined performance criteria are not fulfilled (i.e. nonconformities are present), staff participate in the implementation and recording of corrective actions. The effectiveness of corrective action is monitored.

[H:\Lab\\_Medicine\POCT - Point of Care\ Whole Blood Glucose Meter \(Nova & Roche\)\Trends\Proficiency Testing](H:\Lab_Medicine\POCT - Point of Care\ Whole Blood Glucose Meter (Nova & Roche)\Trends\Proficiency Testing)

- f. Allowable Error:\*may vary depending on sample submitted – see final reports

Provided with the result of the evaluations. Example below

| Specimen | Evaluation Limit | Allowable Limit    |
|----------|------------------|--------------------|
| WBGQ-01  | $\pm 0.7 \%$     | 4.7 – 6.1 mmol/L   |
| WBGQ-02  | $\pm 12.5 \%$    | 13.2 – 17.1 mmol/L |
| WBGQ-03  | $\pm 12.5 \%$    | 17.9 – 23.1 mmol/L |

- i. Problem glucometers will need to be re-verified before return to ward.  
Refer to: [CAPT\\_TDE\\_0300 Validating New/Loaner Glucometer](#)

**REFERENCE:** <https://www.cpsbc.ca/files/pdf/DAP-PT-Manual.pdf>